

Percp/Cyanine5.5 Anti-Mouse Cd25 Monoclonal Antibody For Flow Cytometry - P01590

Item Number: CM0008070

Introduction

High-quality PerCP/Cyanine5.5-conjugated rabbit anti-mouse CD25 monoclonal antibody for flow cytometry. Specifically recognizes mouse CD25 (IL-2 receptor alpha) with cross-reactivity to mouse. Ideal for immune tolerance and Treg cell research.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Mouse CD25 mAb
Remarks/Alias	CD25; IL2r; Ly-43
Species	Mouse
GeneID (Human)	3559
GeneID	16184
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 22-236 of mouse CD25 (NP_032393.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01590
Protein Weight	31kDa
Shipping	Ice bag