

Percp/Cyanine5.5 Anti-Mouse Cd69 Monoclonal Antibody For Flow Cytometry - P37217

Item Number: CM0011666

Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody against mouse CD69. Suitable for flow cytometric detection of activated immune cells. Recognizes a ~22 kDa transmembrane glycoprotein.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Mouse CD69 mAb
Remarks/Alias	AIM; VEA; 5830438K24Rik
Species	Mouse
GeneID (Human)	969
GeneID (Mouse)	12515
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 64-199 of mouse CD69 (NP_001028294.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P37217
Protein Weight	22kDa
Shipping	Ice bag