

Percp/Cyanine5.5 Mouse Anti-Human Cd3 Monoclonal Antibody For Fc - P07766

Item Number: CM0008046



Introduction

PerCP/Cyanine5.5-conjugated mouse monoclonal antibody targeting human CD3 epsilon (CD3E) for flow cytometry. Recognizes the TCR-CD3 complex on T lymphocytes, validated for human samples. Suitable for immunophenotyping and T-cell research applications.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Mouse anti-Human CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
Gene ID (Human)	916
Gene ID	916
Immunogen	Recombinant protein of human CD3.
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag