

Percp/Cyanine5.5 Rabbit Anti-Human Cd27 Monoclonal Antibody For Flow Cytometry - P26842

Item Number: CM0009253



Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody targeting human CD27, a costimulatory receptor on T-cells, NK-cells, and B-cells. Validated for flow cytometry, suitable for immunophenotyping and immune-checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Human CD27 mAb
Remarks/Alias	T14; S152; Tp55; TNFRSF7; S152. LPFS2
Species	Human
Gene ID (Human)	939
Gene ID	939
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-192 of human CD27 (NP_001233.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26842
Protein Weight	29kDa
Shipping	Ice bag