

Percp/Cyanine5.5 Rabbit Anti-Human Cd31/Pecam1 Monoclonal Antibody For Fc - Q08481

Item Number: CM0004654

Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody against human CD31/PECAM1 for flow cytometry. Suitable for detecting CD31 on endothelial cells, platelets, and leukocytes in human samples. SWISS Q08481, molecular weight 83kDa.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Human CD31/PECAM1 mAb
Remarks/Alias	CD31; PECA1; GPIIA'; PECAM-1; endoCAM; CD31/EndoCAM
Species	Human
GeneID (Human)	5175
GeneID	5175
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 28-601 of human CD31/PECAM1 (NP_000433.4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q08481
Protein Weight	83kDa
Shipping	Ice bag