

Percp/Cyanine5.5 Rabbit Anti-Human Cd56/Ncam1 Monoclonal Antibody For Fc - P13591

Item Number: CM0011686



Introduction

Rabbit monoclonal antibody against human CD56/NCAM1, conjugated to PerCP/Cyanine5.5, validated for flow cytometry. Ideal for studying neural cell adhesion, neurite outgrowth, and viral receptor interactions.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Human CD56/NCAM1 mAb
Remarks/Alias	CD56; NCAM; MSK39
Species	Human
Gene ID (Human)	4684
Gene ID	4684
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-413 of human CD56/NCAM1 (NP_851996.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13591
Protein Weight	95 kDa
Shipping	Ice bag