

PerCP/Cyanine5.5 Rabbit Anti-Human CD62L/L-Selectin Monoclonal Antibody For Fc - P14151

Item Number: CM0011411



Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody targeting human CD62L/L-Selectin (P14151). Validated for flow cytometry, this antibody is ideal for studying lymphocyte adhesion and homing mechanisms in human samples.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Human CD62L/L-Selectin mAb
Remarks/Alias	TQ1; LAM1; LEU8; LNHR; LSEL; CD62L; LYAM1; PLNHR; LECAM1
Species	Human
GeneID (Human)	6402
GeneID	6402
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-332 of human CD62L/L-Selectin (NP_000646.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P14151
Protein Weight	42kDa/43kDa
Shipping	Ice bag