

Percp/Cyanine5.5 Rabbit Anti-Human Cd94 Monoclonal Antibody For Flow Cytometry - Q13241

Item Number: CM0010973



Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody targeting human CD94 (KLRD1, NK cell receptor). Designed for flow cytometry to detect human CD94 with specificity. Ideal for immune profiling, NK cell research, and immune checkpoint studies.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Human CD94 mAb
Remarks/Alias	CD94
Species	Human
GeneID (Human)	3824
GeneID	3824
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 32-179 of human CD94. (NP_002253.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13241
Protein Weight	21kDa
Shipping	Ice bag