

Percp/Cyanine5.5 Rabbit Anti-Human/Monkey Cd19 Mab - Cd19 Monoclonal Antibody

Item Number: CM0010835

Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody against human and monkey CD19, validated for flow cytometry. Detects CD19 (63kDa) on B cells. Ideal for immuno-oncology and B-cell research.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Human/Monkey CD19 mAb
Remarks/Alias	CD19
Species	Monkey
Gene ID (Human)	930
Gene ID	102145514
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-292 of monkey CD19.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
Protein Weight	63kDa
Shipping	Ice bag