


Percp/Cyanine5.5 Rabbit Anti-Human/Monkey Cd38 Monoclonal Antibody For Flow Cytometry - P28907

Item Number: CM0008525



Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody targeting human and cynomolgus CD38 (SWISS: P28907). Suitable for flow cytometry, binding the extracellular domain. CD38 is a multifunctional enzyme and receptor involved in calcium signaling and lymphocyte activation.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Human/Monkey CD38 mAb
Remarks/Alias	ADPRC1; cADPR1; ADPRC 1
Species	Human
Gene ID (Human)	952
Gene ID	952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (NP_001766.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P28907
Protein Weight	14kDa/34kDa
Shipping	Ice bag