

Percp/Cyanine5.5 Rabbit Anti-Human/Monkey Cd4 Monoclonal Antibody For Fc - P01730

Item Number: CM0011012



Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody targeting human/monkey CD4, validated for flow cytometry. Ideal for T-helper cell identification, immune profiling, and HIV/HHV-7 research. - P01730

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Human/Monkey CD4 mAb
Remarks/Alias	T4; IMD79; Leu-3; OKT4D; CD4mut
Species	Human
GeneID (Human)	920
GeneID	920
Immunogen	Recombinant protein of human CD4.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01730
Protein Weight	51kDa
Shipping	Ice bag