

Percp/Cyanine5.5 Rabbit Anti-Human/Monkey Cd45 Monoclonal Antibody For Flow Cytometry - P08575

Item Number: CM0012354



Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody targeting human/monkey CD45 (PTPRC, LCA). Validated for flow cytometry, suitable for leukocyte immunophenotyping and T-cell activation studies.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Human/Monkey CD45 mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID (Human)	5788
GeneID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-577 of human CD45 (P08575)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08575
Protein Weight	131kDa/136kDa/141kDa/143kDa/147kDa
Shipping	Ice bag