


Percp/Cyanine5.5 Rabbit Anti-Human/Monkey Cd8A Monoclonal Antibody For Flow Cytometry - P01732

Item Number: CM0009894



Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody against human/monkey CD8a for flow cytometry. Ideal for immune response studies, CTL function analysis, and T-cell subset identification.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Human/Monkey CD8a mAb
Remarks/Alias	CD8; p32; Leu2; CD8alpha
Species	Human
GeneID (Human)	925
GeneID	925
Immunogen	Recombinant fusion protein containing amino acids 22-182 of human CD8A (NP_001759.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus, Rhesus
SWISS	P01732
Protein Weight	21kDa/25kDa/30kDa
Shipping	Ice bag