

Percp/Cyanine5.5 Rabbit Anti-Human/Monkey Pd-1/Cd279 Monoclonal Antibody For Flow Cytometry - Q15116

Item Number: CM0011565



Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody against human and cynomolgus PD-1/CD279 (UniProt Q15116), validated for flow cytometry. Ideal for immune checkpoint research, cancer immunotherapy studies, and T-cell exhaustion analysis.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Human/Monkey PD-1/CD279 mAb
Remarks/Alias	PD1; PD-1; CD279; SLEB2; hPD-1; hPD-I; hSLE1
Species	Human
Gene ID (Human)	5133
Gene ID	5133
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 25-168 of human PD-1/CD279 (NP_005009.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	Q15116
Protein Weight	32 kDa
Shipping	Ice bag