


Percp/Cyanine5.5 Rabbit Anti-Mouse Cd127/Il-7Rα Monoclonal Antibody For Fc - P16872

Item Number: CM0008974



Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody against mouse CD127/IL-7Rα, validated for flow cytometry (FC). Ideal for studying IL-7 receptor signaling and lymphocyte development in immunology research.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Mouse CD127/IL-7Rα mAb
Remarks/Alias	CD127; IL-7Ralpha
Species	Mouse
Gene ID (Human)	3575
Gene ID (Mouse)	16197
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-239 of mouse CD127/IL-7Rα (NP_612564.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P16872
Protein Weight	51kDa
Shipping	Ice bag