

Percp/Cyanine5.5 Rabbit Anti-Mouse Cd19 Monoclonal Antibody For Flow Cytometry (Fc) - P25918

Item Number: CM0012226



Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody targeting mouse CD19, validated for flow cytometry. Ideal for studying B-cell receptor coreceptor function and B-cell differentiation in mouse models.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Mouse CD19 mAb
Species	Mouse
Gene ID (Human)	930
Gene ID	12478
Immunogen	Recombinant protein of mouse CD19.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25918
Protein Weight	60kDa
Shipping	Ice bag