

Percp/Cyanine5.5 Rabbit Anti-Mouse Cd4 Monoclonal Antibody For Flow Cytometry - P06332

Item Number: CM0012175

Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody targeting mouse CD4 (L3T4/Ly-4). Validated for flow cytometry, ideal for immunophenotyping and T-helper cell studies in mouse models.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Mouse CD4 mAb
Remarks/Alias	L3T4; Ly-4
Species	Mouse
Gene ID (Human)	920
Gene ID (Mouse)	12504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-394 of mouse CD4 (NP_038516.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06332
Protein Weight	51kDa
Shipping	Ice bag