

Percp/Cyanine5.5 Rabbit Anti-Mouse Cd44 Monoclonal Antibody For Fc - P15379

Item Number: CM0008772



Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody against mouse CD44 (Ly-24, Pgp-1, HERMES) for flow cytometry (FC). Validated in mouse, suitable for immunophenotyping and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Mouse CD44 mAb
Remarks/Alias	Ly-24; Pgp-1; HERMES
Species	Mouse
GeneID (Human)	960
GeneID	12505
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-176 of mouse CD44 (NP_033981.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P15379
Protein Weight	86kDa
Shipping	Ice bag