

Percp/Cyanine5.5 Rabbit Anti-Mouse Cd68 Monoclonal Antibody For Flow Cytometry - P31996

Item Number: CM0011332

Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody against mouse CD68 (macrosialin, Lamp4, gp110, Scard1). Validated for flow cytometry and intracellular flow cytometry. Suitable for macrophage and dendritic cell research, including phagocytosis and cell trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Mouse CD68 mAb
Remarks/Alias	Lamp4, gp110, Scard1
Species	Mouse
GeneID (human)	968
GeneID	12514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-291 of mouse CD68 (NP_001277987.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, FC (intra)
Cross Reactivity	Mouse
SWISS	P31996
Protein Weight	35kDa
Shipping	Ice bag