

Percp/Cyanine5.5 Rabbit Igg Isotype Control - Monoclonal Antibody

Item Number: CM0011203

Introduction

PerCP/Cyanine5.5-conjugated Rabbit IgG isotype control monoclonal antibody for flow cytometry. Use as a critical negative control to assess nonspecific binding and establish gating strategies in FC experiments. Sourced by CamelBio for IVD raw material requirements.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit IgG isotype control
Immunogen	Synthetic peptide A synthesized peptide derived from rabbit IgG isotype control
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag