

Percp/Cyanine5.5 Rat Anti-Mouse Cd45R/B220 Monoclonal Antibody For Fc - P06800

Item Number: CM0011039

Introduction

PerCP/Cyanine5.5-conjugated rat monoclonal antibody targeting mouse CD45R/B220 for flow cytometry. CD45R/B220, a receptor-type tyrosine-protein phosphatase C, is a key regulator of T-cell activation and lymphocyte development. Ideal for mouse immunophenotyping.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rat anti-Mouse CD45R/B220 mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyl-4
Species	Mouse
GeneID (Human)	5788
GeneID	19264/5788
Immunogen	Recombinant protein of mouse CD45R/B220
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130kDa/140kDa/144kDa
Shipping	Ice bag