

Perilipin-2 Rabbit Mab - Q99541

Item Number: CM0003647

Introduction

High-specificity rabbit monoclonal antibody against human Perilipin-2 (ADFP/ADRP), validated for WB, IHC-P, IF/ICC, and ELISA. Ideal for lipid droplet and metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	Perilipin-2 Rabbit mAb
Remarks/Alias	ADFP; ADRP; Perilipin-2
Species	Human
GeneID	123
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 88-437 of human Perilipin-2 (NP_001113.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q99541
Protein Weight	48kDa
Shipping	Ice bag