

Perilipin-2 Rabbit Monoclonal Antibody For Wb, Elisa - Q99541

Item Number: CM0005161

Introduction

Rabbit monoclonal antibody targeting human Perilipin-2 (ADRP/ADFP), suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for lipid droplet and metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	Perilipin-2 Rabbit mAb
Remarks/Alias	ADFP; ADRP; Perilipin-2
Species	Human
GeneID (Human)	123
GeneID	123
Immunogen	A synthesized peptide derived from human perilipin-2.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99541
Protein Weight	48kDa
Shipping	Ice bag