


Pfdn4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nqp4

Item Number: CM0023277



Introduction

PFDN4 Rabbit Polyclonal Antibody validated in WB, IF/ICC, ELISA; reacts with human, mouse, rat. Targets prefoldin subunit 4 chaperone protein, suitable for protein folding and chaperonin studies.

[Learn More](#)

Parameter	Value
Product Name	PFDN4 Rabbit pAb
Remarks/Alias	C1, PFD4, PFDN4
Species	Human
Gene ID	5203
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-134 of human PFDN4 (NP_002614.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQP4
Protein Weight	15kDa
Shipping	Ice bag