

Pfkfb2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O60825

Item Number: CM0025596



Introduction

Rabbit polyclonal antibody against human PFKFB2, validated for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for studying glycolysis regulation in heart tissue.

[Learn More](#)

Parameter	Value
Product Name	PFKFB2 Rabbit pAb
Remarks/Alias	PFK-2/FBPase-2; PFKFB2
Species	Human
GeneID (Human)	5208
GeneID	5208
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 446-505 of human PFKFB2 (NP_006203.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60825
Protein Weight	58kDa
Shipping	Ice bag