


Pfkfb3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q16875

Item Number: CM0006089



Introduction

High-quality PFKFB3 Rabbit Monoclonal Antibody validated for WB, IHC-P, and ELISA. Reacts with human and mouse PFKFB3 (60kDa), a key regulator of glycolysis. Ideal for cancer metabolism research.

[Learn More](#)

Parameter	Value
Product Name	PFKFB3 Rabbit mAb
Remarks/Alias	PFK2, IPFK2, iPFK-2, PFKFB3
Species	Human
Gene ID (Human)	5209
Gene ID	5209
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16875
Protein Weight	60kDa
Shipping	Ice bag