

Pgk1 Rabbit Monoclonal Antibody For Wb, Elisa - P00558

Item Number: CM0002571



Introduction

High-quality rabbit monoclonal antibody targeting human PGK1 (Phosphoglycerate kinase 1), suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for glycolysis, mitochondrial function, and sperm motility studies.

[Learn More](#)

Parameter	Value
Product Name	PGK1 Rabbit mAb
Remarks/Alias	PGKA; MIG10; HEL-S-68p; PGK1
Species	Human
Gene ID (Human)	5230
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00558
Protein Weight	45kDa
Shipping	Ice bag