

Pgm1 Rabbit Pab For Wb, Elisa - P36871

Item Number: CM0028413



Introduction

Rabbit polyclonal antibody targeting human Phosphoglucomutase-1 (PGM1). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for glucose metabolism and glycogen storage disease research.

[Learn More](#)

Parameter	Value
Product Name	PGM1 Rabbit pAb
Remarks/Alias	CDG1T; GSD14; PGM1
Species	Human
Gene ID (Human)	5236
Gene ID	5236
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PGM1 (NP_002624.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36871
Protein Weight	61kDa
Shipping	Ice bag