

Phd3 Rabbit Mab For Wb, Ihc-P, If/Icc, Elisa - Q9H6Z9

Item Number: CM0003927

Introduction

High-quality PHD3 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, ELISA; reacts with human and mouse. Essential tool for hypoxia, cancer, and ferroptosis research.

[Learn More](#)

Parameter	Value
Product Name	PHD3 Rabbit mAb
Remarks/Alias	PHD3; HIFPH3; HIFP4H3
Species	Human
GeneID (Human)	112399
GeneID	112399
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PHD3 (NP_071356.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H6Z9
Protein Weight	27kDa
Shipping	Ice bag