

Phf21A Rabbit Polyclonal Antibody For Wb And Elisa - Q96Bd5

Item Number: CM0025241



Introduction

Validated PHF21A rabbit polyclonal antibody for Western blot and ELISA. Detects human, mouse, rat PHF21A, a BHC corepressor component involved in neuronal gene repression and chromatin modification.

[Learn More](#)

Parameter	Value
Product Name	PHF21A Rabbit pAb
Remarks/Alias	BHC80; NEDMS; BM-006; IDDBCS; PHF21A
Species	Human
Gene ID	51317
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 461-680 of human PHF21A (NP_001095272.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96BD5
Protein Weight	75kDa
Shipping	Ice bag