


Phospho-A-Synuclein-S129 Rabbit Polyclonal Antibody For Wb, Elisa - P37840

Item Number: CM0001863



Introduction

Phospho- α -Synuclein-S129 Rabbit Polyclonal Antibody validated for WB and ELISA in human, mouse, rat. Targets phosphorylated serine 129 of alpha-synuclein (SNCA), key for Parkinson's disease research. - P37840

[Learn More](#)

Parameter	Value
Product Name	Phospho- α -Synuclein-S129 Rabbit pAb
Remarks/Alias	PD1; NACP; PARK1; PARK4; Phospho- α -Synuclein-S129
Species	Human
Gene ID (Human)	6622
Gene ID	6622
Immunogen	A synthetic phosphorylated peptide around S129 of human alpha-Synuclein (NP_000336.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37840
Protein Weight	14kDa
Shipping	Ice bag