


Phospho-A-Synuclein-Y133 Rabbit Pab For Wb, Elisa - P37840

Item Number: CM0000937



Introduction

Phospho- α -Synuclein-Y133 Rabbit polyclonal antibody validated for WB and ELISA. Reacts with human and mouse samples. Useful for investigating synaptic activity, dopamine neurotransmission, and neurodegenerative disorders like Parkinson's and Alzheimer's disease.

[Learn More](#)

Parameter	Value
Product Name	Phospho- α -Synuclein-Y133 Rabbit pAb
Remarks/Alias	PD1; NACP; PARK1; PARK4; Phospho- α -Synuclein-Y133
Species	Human
GeneID (human)	6622
GeneID	6622
Immunogen	Synthetic phospho-specific peptide corresponding to residues surrounding Y133 of human α -Synuclein
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P37840
Protein Weight	14kDa
Shipping	Ice bag