


Phospho-Acc1 (Ser79) Rabbit Polyclonal Antibody For Wb, Elisa - Q13085

Item Number: CM0001118



Introduction

Rabbit polyclonal antibody targeting phospho-ACC1 (Ser79) validated for WB and ELISA; reacts with human, mouse, and rat ACC1; suitable for fatty acid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ACC1-S79 Rabbit pAb
Remarks/Alias	ACC; ACAC; ACC1; ACCA; Acac1; hACC1; ACACAD; ACCalpha; ACACalpha
Species	Human
Gene ID	31
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S79 of human ACC1 (NP_942133.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13085
Protein Weight	266kDa
Shipping	Ice bag