


Phospho-Acc1/Acc2-S79 Rabbit Polymab® For Wb, If/Icc, Elisa - Q13085 / 000763

Item Number: CM0000751



Introduction

Phospho-specific rabbit recombinant monoclonal antibody detecting acetyl-CoA carboxylase 1/2 phosphorylated at Ser79. Validated for WB, IF/ICC, ELISA in human, mouse, rat. Ideal for fatty acid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ACC1/ACC2-S79 Rabbit PolymAb®
Remarks/Alias	ACC; ACAC; ACC1; ACCA; Acac1; hACC1; ACACAD; ACCalpha; ACACalpha; Phospho-Acetyl CoA Carboxylase-S79
Species	Human
Gene ID (Human)	31/32
Gene ID	31/32
Immunogen	Synthetic phosphorylated peptide around S79 of human Acetyl CoA Carboxylase.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13085/000763
Protein Weight	265kDa
Shipping	Ice bag