


Phospho-Acly-S455 Rabbit Mab For Wb Ihc-P Elisa - P53396

Item Number: CM0001626



Introduction

Phospho-specific rabbit monoclonal antibody targeting human ACLY Ser455, validated for Western blot (WB), immunohistochemistry (IHC-P), and ELISA. Cross-reacts with mouse and rat. Ideal for lipid metabolism and signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ACLY-S455 Rabbit mAb
Remarks/Alias	ACL; ATPCL; CLATP; Phospho-ACLY-S455
Species	Human
GeneID (Human)	47
GeneID	47
Immunogen	A synthetic phosphorylated peptide around S455 of human ACLY (NP_001087.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53396
Protein Weight	121kDa
Shipping	Ice bag