


Phospho-Acly-S455 Rabbit Pab For Wb, Ihc-P, Ip, Elisa - P53396

Item Number: CM0000515



Introduction

High-quality rabbit polyclonal antibody against human ACLY phosphorylated at Ser455. Validated for WB, IHC-P, IP, and ELISA. Cross-reacts with mouse and rat. Ideal for lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ACLY-S455 Rabbit pAb
Remarks/Alias	ACL; ATPCL; CLATP; Phospho-ACLY-S455
Species	Human
Gene ID (Human)	47
Gene ID	47
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S455 of human ACLY (NP_001087.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53396
Protein Weight	121kDa
Shipping	Ice bag