

Phospho-Akt (S473) Rabbit Monoclonal Antibody For Wb, Elisa - P31749 / P31751 / Q9Y243

Item Number: CM0001587



Introduction

Phospho-Akt (S473) rabbit monoclonal antibody validated for WB and ELISA. Specifically detects endogenous levels of Akt only when phosphorylated at Ser473 of human, mouse, and rat origin. Ideal for studying PI3K/AKT signaling in cancer and metabolism.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Akt-S473 Rabbit mAb
Remarks/Alias	AKT1/AKT2/AKT3; Phospho-Akt-S473
Species	Human
GeneID (Human)	207/208/10000
GeneID	207/208/10000
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S473 of human AKT1 (NP_005154.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749/P31751/Q9Y243
Protein Weight	48kDa/55kDa/51kDa/54kDa
Shipping	Ice bag