

Phospho-Akt-S473 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P31749 / P31751 / Q9Y243

Item Number: CM0000663



Introduction

Phospho-specific rabbit polyclonal antibody recognizing Akt1, Akt2, and Akt3 phosphorylated at Ser473. Validated for Western blot, immunohistochemistry, and ELISA in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Akt-S473 Rabbit pAb
Remarks/Alias	AKT1/AKT2/AKT3; Phospho-Akt-S473
Species	Human
Gene ID (Human)	207/208/10000
Gene ID	207/208/10000
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749/P31751/Q9Y243
Protein Weight	48kDa/55kDa/51kDa/54kDa
Shipping	Ice bag