

Phospho-Akt-Y315 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P31749 / P31751 / Q9Y243

Item Number: CM0001715



Introduction

Phospho-Akt-Y315 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat Akt phosphorylated at Y315. Suitable for studying the PI3K/Akt pathway in cancer and metabolism.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Akt-Y315 Rabbit pAb
Remarks/Alias	AKT1/AKT2/AKT3; Phospho-Akt-Y315
Species	Human
Gene ID (Human)	207/208/10000
Gene ID	207/208/10000
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around Y315 of human AKT1 (NP_001014431.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749/P31751/Q9Y243
Protein Weight	48kDa/55kDa/51kDa/54kDa
Shipping	Ice bag