

Phospho-Akt1-S124 Rabbit Monoclonal Antibody For Wb And Elisa - P31749

Item Number: CM0001936



Introduction

Phospho-specific rabbit monoclonal antibody targeting human AKT1 phosphorylated at serine 124. Validated for Western blot and ELISA. Suitable for studying AKT signaling in cancer, metabolism, and cell survival pathways.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-AKT1-S124 Rabbit mAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; Phospho-AKT1-S124
Species	Human
Gene ID (Human)	207
Gene ID	207
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S124 of human Akt1 (P31749).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P31749
Protein Weight	56kDa
Shipping	Ice bag