

Phospho-Akt1-S129 Rabbit Mab For Wb And Elisa - P31749

Item Number: CM0001918



Introduction

Rabbit monoclonal antibody against AKT1 phosphorylated at Ser129, validated for western blot and ELISA with human samples. Ideal for studying AKT signaling in cancer, metabolism, and cell survival research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AKT1-S129 Rabbit mAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; Phospho-AKT1-S129
Species	Human
Gene ID	207
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S129 of human Akt1 (P31749)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P31749
Protein Weight	56kDa
Shipping	Ice bag