

Phospho-Akt1-S473 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P31749

Item Number: CM0001804



Introduction

Phospho-specific rabbit polyclonal antibody against human AKT1 phosphorylated at S473. Suitable for WB, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for studying AKT signaling in cancer and metabolism.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AKT1-S473 Rabbit pAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; Phospho-AKT1-S473
Species	Human
Gene ID (Human)	207
Gene ID	207
Immunogen	Synthetic phosphorylated peptide around S473 of human Akt1 (NP_005154.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749
Protein Weight	56kDa
Shipping	Ice bag