

Phospho-Akt1-T308 Rabbit Pab - P31749

Item Number: CM0001290



Introduction

Phospho-AKT1-T308 rabbit polyclonal antibody validated for WB and ELISA. Detects phosphorylated human AKT1 (Thr308) with cross-reactivity to mouse and rat. Ideal for studying insulin signaling, cancer biology, and cell survival pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AKT1-T308 Rabbit pAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; Phospho-AKT1-T308
Species	Human
GeneID (Human)	207
GeneID	207
Immunogen	Synthetic peptide phosphorylated around T308 of human Akt1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749
Protein Weight	56kDa
Shipping	Ice bag