

Phospho-Akt1-T308+Akt2-T309+Akt3-T305 Rabbit Polyclonal Antibody For Wb, Elisa - P31749 / P31751 / Q9Y243

Item Number: CM0001750



Introduction

Phospho-AKT1-T308+AKT2-T309+AKT3-T305 rabbit polyclonal antibody validated for Western blot and ELISA. Detects human and mouse phosphorylated AKT isoforms, suitable for cell signaling and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AKT1-T308+AKT2-T309+AKT3-T305 Rabbit pAb
Remarks/Alias	AKT1/AKT2/AKT3; Phospho-AKT1-T308+AKT2-T309+AKT3-T305
Species	Human
GeneID (Human)	207/208/10000
GeneID	207/208/10000
Immunogen	Synthetic phosphorylated peptide around T308 of human AKT1 (NP_005154.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P31749/P31751/Q9Y243
Protein Weight	48kDa/55kDa/51kDa/54kDa
Shipping	Ice bag