

Phospho-Akt1-T450 + Akt2-T451 + Akt3-T447 Rabbit Mab For Wb, Ihc-P, Elisa - P31749 / P31751 / Q9Y243

Item Number: CM0000773



Introduction

Phospho-specific rabbit monoclonal antibody detecting AKT1 phosphorylated at Thr-450, and AKT2/AKT3 at equivalent sites. Suitable for WB, IHC-P, and ELISA; cross-reacts with human, mouse, and rat AKT isoforms. Ideal for studying PI3K/AKT signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AKT1-T450 + AKT2-T451 + AKT3-T447 Rabbit mAb
Species	Human
Gene ID (Human)	207/208/10000
Gene ID	207/208/10000
Immunogen	Synthetic peptide (A synthetic phosphorylated peptide around T450 of human Akt)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749/P31751/Q9Y243
Protein Weight	48KDa/51KDa/55KDa
Shipping	Ice bag