


Phospho-Akt1S1/Pras40-T246 Rabbit Monoclonal Antibody For Wb, Elisa - Q96B36

Item Number: CM0000766



Introduction

Phospho-specific rabbit monoclonal antibody targeting human AKT1S1/PRAS40 phosphorylated at Thr246. Validated for WB and ELISA. Ideal for studying mTORC1 signaling, cancer metabolism, and insulin-responsive pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AKT1S1/PRAS40-T246 Rabbit mAb
Remarks/Alias	Lobe; PRAS40
Species	Human
GeneID (Human)	84335
GeneID	84335
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T246 of human AKT1S1/PRAS40 (NP_001265088.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96B36
Protein Weight	27 kDa
Shipping	Ice bag