


Phospho-Akt1S1/Pras40-T246 Rabbit Pab For Wb, Elisa - Q96B36

Item Number: CM0001273



Introduction

This phospho-specific rabbit polyclonal antibody detects AKT1S1/PRAS40 phosphorylated at T246 in human samples via WB and ELISA. Ideal for mTORC1 signaling and insulin/nutrient response studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AKT1S1/PRAS40-T246 Rabbit pAb
Remarks/Alias	Lobe; PRAS40; Phospho-AKT1S1/PRAS40-T246
Species	Human
Gene ID (Human)	84335
Gene ID	84335
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96B36
Protein Weight	27kDa
Shipping	Ice bag