


Phospho-Akt2 Specific-S474 Rabbit Pab - P31751

Item Number: CM0001731



Introduction

Phospho-specific rabbit polyclonal antibody detecting AKT2 phosphorylated at Ser474, validated for WB and ELISA. Ideal for investigating insulin signaling, glucose metabolism, and cancer cell survival pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AKT2 specific-S474 Rabbit pAb
Remarks/Alias	PKBB; PRKBB; HIHGH; PKBBETA; RAC-BETA; Phospho-AKT2 specific-S474
Species	Human
GeneID (human)	208
GeneID	208
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P31751
Protein Weight	56kDa
Shipping	Ice bag