

Phospho-ALK-Y1096 Rabbit Pab For Wb, Elisa - Q9Um73

Item Number: CM0001181



Introduction

Phospho-ALK-Y1096 Rabbit pAb: Detects phosphorylated ALK at tyrosine 1096 in human samples. Validated for WB and ELISA. Key for studying ALK signaling in neurogenesis, energy balance, and oncogenic pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ALK-Y1096 Rabbit pAb
Remarks/Alias	ALK1; CD246; NBLST3; Phospho-ALK-Y1096
Species	Human
GeneID (Human)	238
Immunogen	Synthetic peptide; A synthetic phosphorylated peptide around Y1096 of human ALK (NP_004295.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UM73
Protein Weight	176kDa
Shipping	Ice bag