

Phospho-Alt-Rpl36-S140 Rabbit Polyclonal Antibody - Alt-Rpl36

Item Number: CM0001784

Introduction

Phospho-alt-RPL36 (S140) rabbit polyclonal antibody validated for Western blot, dot blot, and ELISA in human samples. Ideal for detecting phosphorylation events in alt-RPL36 to support signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-alt-RPL36-S140 Rabbit pAb
Species	Human
Immunogen	Synthetic Phosphorylated Peptide (A synthetic phosphorylated peptide around S140 of alt-RPL36)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human
Shipping	Ice bag